

NUWEIGH®

LOG-984 WHEEL CHAIR SCALE



Carry handles

4x High quality load cells

Beam weight of 8.1kg

Weighs in kg & lb

Entry and exit ramps

Lightweight & Portable

INDICATOR

PRE-TARE FUNCTION

Allows the user to manually enter the wheel chair weight prior to weighing the patient.

PRE TARE

BMI

TARE

HOLD

UNIT



Model

Capacity

Division

Part No.


LOG-984

300kg

100g

#1712

DIMENSIONS



ALWAYS FOLLOW THE DIRECTIONS FOR USE