

NUWEIGH®

XAM710 Professional Personal Wireless Scale

NUWEIGH's Professional Personal Scale uses four precision load cells for high accuracy readings. The wireless display makes the XAM710 a perfect scale for medical practitioners

Non-slip platform

Reads in both kgs or lbs

9V power adapter or
4 x A4 battery operation

Hold and BMI function

Auto OFF function for
power saving

Lightweight and
portable structure

Four precision load cells



Model	Capacity	Resolution	Platform Size	Inv. Number
XAM710	250kg	100g	350 x 345mm	36608

Ideal for medical practitioners

Wireless 26mm display can be located up to 30 metres from the scale platform



Wireless Display

Connect up to three scale bases to one display!

