NUWEIGH

CHAIR SCALE CHR-691

FEATURES

- Adjustable arm rests allow for easy entry/exit from the chair
- 300kg capacity
- 100gram accuracy
- Strong fiberglass chair
- BMI function
- Tubular design with no sharp edges
- Lift up swivel footrests for comfort and easier patient transfer
- Level indicator for consistent accuracy
- Large 26mm LCD display
- · Built in rechargeable battery
- All cabling is protected
- Ergonomic 930mm handle height
- · 4 Swivel castor wheels with locks
- · Hygienic ABS seat
- Scale approved by Therapeutic Goods Administration





#1533



300kg x 100g

#200237