

# NUWEIGH®



## MS3830 Wheel Chair Scale

### Features

NUWEIGH's wheelchair scale is fitted with a robust handle and wheels for easy transportation and positioning. The integrated ramps allow easy accessibility when entering and exiting the platform

- Large 25mm LCD display
- Pre-tare function for easy wheel chair weighing
- Adjustable feet
- 1.5v AAA battery or 240V
- BMI function
- 300kg capacity x 100g accuracy
- TGA Approved



Code	Model	Capacity	Display	Dimensions
#4673	MS3830	300kg x 100g	25mm LCD screen	740mm x 1140mm



Optional left and right stability rails/  
supporting brace



Integrated handle and  
wheels for easy movement