

# NUWEIGH®

## NUWEIGH® LOG842 Electronic Medical Scale


*Low weight, extremely robust construction and easy-to-use display.*



### FEATURES

- › 300kg Capacity
- › Four (4) load cell design
- › Large platform size of 310mm x 300mm
- › Large 26mm LCD design
- › Lightweight and portable
- › Auto power off
- › Weighs in both kg/lb
- › 9V Battery or 240V operation
- › Power adapter included

### FUNCTIONS



### DIMENSIONS



#### 300kg Capacity

Four highly accurate and reliable load cells



#### Lightweight Structure

Robust, yet portable device is ideal for both personal use and medical professionals



*\*On Request*

#### Calibration Certificate

Professionally calibrated by a trained scale technician



*\*Optional Extra*

#### Carry Bag

Purpose designed carry bag protects and increases the life span of the LOG842 Scale

**NWS**  
NUWEIGH

**NUWEIGH® AUSTRALIA**

NATIONAL SERVICE AND SALES NSW | VIC | ACT | QLD | NT | SA | WA | TAS

1300 885 746 | sales@nuweigh.com.au | www.nuweigh.com.au

Newcastle Weighing Services EST.1976 | ABN 36 001 344 430