

# NUWEIGH

## PROFESSIONAL PERSONAL SCALE

### XAM420

NUWEIGH'S Professional Personal Scale uses four precision load cells for high accuracy readings.



**ACCURATE  
WEIGHT  
READINGS**

#### FUNCTIONS

HOLD

TARE

BMI

LCD



Dual sided display at 26mm allows both patient and practitioner to view reading

#### "2 scales in 1" Easily weigh mother and infant seperately


This scale will weigh adult 250kg by 100g then reset to weigh an infant 20kg by 10g whilst held by adult.




#### FEATURES

- Four precision loadcells
- Auto OFF function for power saving
- Hold and BMI function
- Non-slip platform
- Reads in both kgs or lbs
- 9V power adaptor or 4 x A4 battery
- Lightweight and portable structure

#### DIMENSIONS



Extra large weighing platform



Capacity	Divisions	Part No.
250kg	100g	#36608