


NUWEIGH

Stainless Steel Test Masses



Pictured: 5kg/10kg/20kg rectangular weights



FEATURES

- E2 class polished stainless steel
- F1 class polished stainless steel
- M1 class polished stainless steel
- 1mg to 20kg weights available
- Adheres to OIML R111 standard
- NATA and REG13 Certificates available
- Manufacturing under strict quality control
- High quality design

NWS
NUWEIGH

CAPACITIES AVAILABLE

1mg	1g
2mg	2g
5mg	5g
10mg	10g
20mg	20g
50mg	50g
100mg	100g
200mg	200g
500mg	500g
1kg	2kg
5kg	10kg
20kg	

Heavy duty aluminium boxes available for individual weights or sets



NUWEIGH AUSTRALIA • 1300 885 746 • www.nuweigh.com.au

NWS PO Box 326 Wickham NSW 2293 • Established in 1976 • ABN 36 001 344 430

NWS
NUWEIGH