


NUWEIGH CHARGER

HM200PW

NWS
NUWEIGH



The height of person plays an important role in treatment, especially Body Mass Index (BMI). BMI is a calculated value which allows you to define how well the mass of the human body corresponds to the height and consequently a degree of over or under weight.

Using BMI you can try to estimate the risk of developing an illness due to excess or low weight. Easy to assemble and simple to use HM200PW makes measuring height an easy routine job.

Product Description

HM200PW is a height measuring rod for mounting on the wall.

- Measuring range – 14 ~ 205cm
- Division – 1 mm
- Manual type



Product Specifications #25326

Measuring range	14 ~ 205 cm
Graduation	1 mm
Dimensions	350 x 400 x 2430 mm
Weight	3.0 kg
Power supply	–
Features	Self-assembly height rod to be fixed on the wall

ALWAYS FOLLOW THE DIRECTIONS FOR USE

NUWEIGH AUSTRALIA • 1300 885 746 • www.nuweigh.com.au

NWS PO Box 326 Wickham NSW 2293 • Established in 1976 • ABN 36 001 344 430

NWS
NUWEIGH