

NUWEIGH®

LOG244

BABY SCALE

2.6kg Lightweight design for easy portability

Ultra slim profile

4 Load cell design

9V Power adaptor or 240V battery operation



Light and robust construction

Weigh in both lb/kg

Removable head-rests

FUNCTIONS

- ▶ Auto Hold
- ▶ Auto Off
- ▶ Zero
- ▶ Tare

Baby movement compensation technology

Large, easy to read LCD display

ALWAYS FOLLOW THE DIRECTIONS FOR USE

INCLUDES QUALITY CARRY CASE



EASY TO READ 26mm DISPLAY




OPTIONAL ADJUSTABLE HEIGHT ROD

MEASURING RANGE: 35cm - 80cm, 1mm increments



DIMENSIONS



CAPACITY	RESOLUTION	PART NO.
20kg	10g	#1709



NUWEIGH® AUSTRALIA

NATIONAL SERVICE AND SALES NSW | VIC | ACT | QLD | NT | SA | WA | TAS

1300 885 746 | sales@nuweigh.com.au | www.nuweigh.com.au

Newcastle Weighing Services EST.1976 | ABN 36 001 344 430