

NUWEIGH®

LOG344 Baby Scale



**Accurate measurements
with baby movement
compensation technology**

Features

- Deep sided baby cradle
- Cradle can be easily detached for cleaning
- Dual range accuracy
- Hold function
- Large 26mm LCD display
- Remote indicator
- 6 x AAA batteries and 240V adapter



Code	Model	Dimensions	Weight	Capacity	Increments
#37253	LOG344	290mm(L) x 550mm(W)	3.70kg	6kg/15kg	2g/5g

Weigh your baby with the precision of the LOG344 Baby Scale. Specifically designed to weigh undersized or premature babies.


OPTIONAL EXTRA - NUWEIGH Height Rod



Complete your compact baby weighing station by adding a NUWEIGH Height Rod attachment.

Measuring range: 35cm - 80cm
in 1mm divisions

A NUWEIGH Height Rod is designed to attach seamlessly to your LOG344



LOG344 BABY SCALE

ALWAYS FOLLOW THE DIRECTIONS FOR USE

National Sales and Service NSW | VIC | ACT | QLD | NT | SA | WA | TAS
1300 885 746 | www.nuweigh.com.au

NWS
NUWEIGH