

NUWEIGH®

LOG507 PROFESSIONAL MEDICAL SCALE

AAA Battery and
240V Operation

NUWEIGH®
Trademark

Light & Robust
Construction

FUNCTIONS

- ▶ Auto Off
- ▶ Tare
- ▶ Zero
- ▶ Hold

BMI
Function

ERGONOMIC CARRY HANDLE



LARGE ADJUSTABLE FEET



LARGE 25mm LCD DISPLAY



NON-SLIP RUBBER SURFACE



OPTIONAL CARRY BAG




INTEGRATED CARRY HANDLE SUITS MOBILE
MEDICAL PROFESSIONALS



The easy-to-read LCD display and precision weight readings make the **LOG507 PROFESSIONAL MEDICAL SCALE** the perfect weighing solution for all medical facilities and health care workers. Weighing bariatric patients is safe and easy, thanks to the non-slip platform, die-cast housing and 300kg capacity.

DIMENSIONS



**HIGH
CAPACITY**



**CARRY
HANDLE**



**CERTIFICATE
AVAILABLE**

ALWAYS FOLLOW THE DIRECTIONS FOR USE

NWS
NUWEIGH

NUWEIGH® AUSTRALIA

NATIONAL SERVICE AND SALES NSW | VIC | ACT | QLD | NT | SA | WA | TAS

1300 885 746 | sales@nuweigh.com.au | www.nuweigh.com.au

Newcastle Weighing Services EST.1976 | ABN 36 001 344 430