

NUWEIGH®

LOG672



Features:

Design

Sturdy, robust design ensuring safety and longevity

BMI Function

Simple to use BMI function with scale and height rod

Safety Weigh Rail

Safety rail ensures frail or barometric patients can balance themselves while being weighed.

Battery or Powerpack

5 x 1.5v Battery or 240v Powerpack


RS232 Connection

Pre-Tare Function

Base Size

Large base size
550mm x 550mm

Numerical Keyboard



300kg

300kg weight capacity

Portable

Easily transportable with coaster wheels

Code	Model	Capacity	Resolution
#4154	LOG672	300kg	100g

ALWAYS FOLLOW THE DIRECTIONS FOR USE

NWS
NUWEIGH

NATIONAL SALES AND SERVICE NSW | VIC | ACT | QLD | NT | SA | WA | TAS
1300 885 746 | sales@nuweigh.com.au | www.nuweigh.com.au | EST.1976 | ABN 36 001 344 430