

# NUWEIGH®

## LOG842 Electronic Medical Scale

Low weight, extremely robust construction and easy-to-use display.



### Features

- › 300kg Capacity
- › Four (4) load cell design
- › Large platform size of 310mm x 300mm
- › Large 26mm LCD design
- › Lightweight and portable
- › Auto power off
- › Weighs in both kg/lb
- › 9V Battery or 240V operation
- › Power adapter included

### Functions



### Dimensions



### 300kg Capacity

Four highly accurate and reliable load cells



### Lightweight Structure

Robust, yet portable device is ideal for both personal use and medical professionals



*\*On Request*

### Calibration Certificate

Professionally calibrated by a trained scale technician



*\*Optional Extra*

### Carry Bag

Purpose designed carry bag protects and increases the life span of the LOG842 Scale

**ALWAYS FOLLOW THE DIRECTIONS FOR USE**