

NUWEIGH®

LOG909 Scale and Height Rod

Features

- Heavy-duty rugged aluminium design
- Convenient slide out height rod in pole
- Easy to use BMI function with graph
- Wheels for portability
- Nonslip rubber coated platter
- Capacity 250kg
- Accuracy 100g
- Rechargeable battery or 240V operation
- All NUWEIGH scales are calibrated by expert weighing technicians with precision masses



Purpose built for heavy-duty use in hospitals and medical environments



Code	Model	Capacity	Accuracy	Dimensions
#4726	LOG909	250kg	100g	Platter: 370mm x 370mm Height: 1020mm

ALWAYS FOLLOW THE DIRECTIONS FOR USE



Measure height, weight and BMI with the **LOG909**