

NUWEIGH®

STAINLESS STEEL TEST MASSES



CLASSES

- E2 class polished stainless steel
- F1 class polished stainless steel
- M1 class polished stainless steel

FEATURES

- 1mg to 20kg weights available
- Adheres to OIML R111 standard
- Manufactured under strict quality control
- High quality design

CAPACITIES AVAILABLE

1mg	1g
2mg	2g
5mg	5g
10mg	10g
20mg	20g
50mg	50g
100mg	100g
200mg	200g
500mg	500g
1kg	2kg
5kg	10kg
20kg	

Heavy duty aluminium boxes available for individual weights or sets



**NATA
& REG13
REPORTS
AVAILABLE**

***CAN MANUFACTURE CUSTOM
SIZE AND CAPACITY**